Precautions

1. Clean exposed parts with a soft, slightly damp cloth. To remove stains, use a mild soap. Never use detergents, excess water, treated cloths, harsh cleaning agents, or sprays. Do not attempt to lubricate or open scale casing.
2. Although your scale is designed to be maintenance-free, it contains sensitive electronic parts. Avoid rough treatment. Do not drop, kick, or jump on it. Treat it with care to ensure the best performance.
3. This scale is designed to read weights between 4.4 and 400 lbs and 2 and 180 kgs. To prevent damage, do not place more than 400 lbs / 180 kgs on it.
4. Do not store scale where you store cleaning chemicals. The vapors from some household products may affect the electronic components of your scale. Do not store scale on its side.
5. This scale is an extremely sensitive weighing device. In order to prevent run down of the battery, do not store anything on the scale.
6. Do not dispose of batteries in fire. Batteries may explode or leak. Remove the battery if the scale will not be used for a long period of time.

Caution: Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Taylor Precision Products, Inc.
2220 Entrada Del Sol, Suite A
Las Cruces, New Mexico 88001 USA

©2014 Taylor Precision Products, Inc. and its affiliated companies, all rights reserved. Taylor® and Leading the Way in Accuracy® are a registered trademarks of Taylor Precision Products, Inc. and its affiliated companies. All rights reserved.

Made to our exact specifications in China.
Not legal for trade.
www.taylorusa.com

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Leading the Way in Accuracy  

Instruction Manual

Thank you for purchasing a Taylor   Lithium Electronic scale. All of our precision instruments for weight measurement are engineered and designed to meet the highest quality standards ...to assure you years of uncompromising accuracy and consistently dependable, convenient performance.

These instructions contain important and useful information about this unit's operation. Please read instructions thoroughly to ensure that you have the full benefit of all the unit's features. Keep these instructions handy for future reference.

General Operation

1. There is a strip under the Lithium battery cover to prevent battery drain in shipping. Please remove the strip before continuing with these instructions (you may need to remove the battery compartment cover to avoid breaking the strip). Also, some models have a static cling label on the lens to prevent scratching. Please remove before use.
2. Your scale has been set at the factory to read weight in pounds. You may select kilograms units of measurement by moving the slide switch on the bottom of the scale.
3. Place scale on flat, hard surface. Carpeted or uneven floors may affect accuracy.
4. Step on the scale (disregarding the reading) and step off. This process insures the scale is set and is now ready for use. At all other times step straight onto the scale.
5. Step onto the scale platform and remain still. The display will count up to your weight. Your weight is now displayed.
6. The scale will turn off automatically.

Battery Information

Your Taylor   scale has a lithium battery installed to provide extra long energy life for your scale. Under heavy usage (Health Club, Wellness Clinic, etc.) the battery may eventually need to be replaced.

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Battery Replacement

- Use 1- 3 volt CR2032 lithium battery (as installed).
- Remove screw (or screws) from the battery compartment cover on the bottom of the scale and remove the cover.
- Remove the old battery from the compartment.
- Place the new battery with the “+” sign up into the battery compartment.
- Replace the batteries compartment cover and screw(s).

NOTE: Please recycle or dispose of batteries per local regulations.

WARNING: Batteries may pose a choking hazard. As with all small items, do not let children handle batteries. If swallowed, seek medical attention immediately.

PRECAUTION: Do not dispose of batteries in fire. Batteries may explode or leak. Remove the battery if the scale will not be used for a long period of time.

Low Battery

Your Taylor® scale is equipped with a low battery indicator. It will display “Lo” or readings will grow dim or become irregular if the battery needs replacement.

Trouble Shooting

In Case of Difficulty...

- Place your feet in the same position every time you weigh yourself. Stand still during measurement.
- Weigh yourself without clothing whenever possible. Clothes can vary in weight.
- Move the scale to a location where it won't be bumped. If your scale is bumped while being activated or while in use, an inaccurate reading may occur.
- For the most consistent readings, weigh yourself at the same time of day. Allow for weight fluctuations due to food or drink in your system.
- “Err0” Initialization error. To correct, press on the platform. The scale will turn on and off. Repeat measurement.
- “Err1” Instability Error. Ensure the scale is on a flat, stable surface. Step off and back onto the scale, standing still while your weight computes.
- The maximum capacity of this scale is 400 lbs / 180 kgs. Should “O-Ld” appear on the display during a weighing, this indicates that the maximum weighing capacity of the scale has been exceeded. Remove the weight immediately; otherwise, permanent damage may occur.

Reaching Target Weight

Maintaining proper weight is essential to healthy living. Obesity has been linked to a higher risk of heart disease, diabetes and some forms of cancer. The Body Mass Index (BMI) is the most widely accepted measurement of weight according to health professionals.

Use the chart below to determine your own BMI by matching your height in the left hand column with your weight in the center. Anyone, male or female, with a BMI of 25 to 29 is considered over weight and a body mass index of 30 is considered obese. Remember to eat right and exercise.

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